



A COMPLETE VETERINARY SPECIFIC DIGITAL IMAGING SYSTEM





DelWorks DR

DelWorks Veterinary DR is a complete digital package specifically designed for the veterinary imaging. High resolution images are captured with the 17in x 17in TOSHIBA FDX4343R, Cesium Iodide (CsI) flat panel detector.

The user-friendly DelWorks Veterinary specific acquisition software is easy to operate and features a wide range of tools for a variety of image manipulations. Images display in less than 1 second after exposure on the TFT color LCD clinical review display with touch screen interface. This DICOM compliant, high definition display provides more accurate gray scales and better contrast than standard monitors.

FDX4343R

Automatic image optimization is achieved with the integrated *Advanced Image Processing* software, so you can get the image you need fast to make your diagnosis quickly and efficiently.





···· 12

DelWorks Veterinary Image Acquisition & Processing Software



Simple to use software that guides you every step of the way.

Begin by planning your exam. Choose the animal and area(s) of interest with DelWorks anatomical selection tool. Pre-programmed plan studies included for dogs, cats, horses and several small pets.



The unique Animal Positioning Guide ensures accuracy with step by step instructions and visual aids. The generator panel displays all values and settings recommended for a specific examination.



View your image on the preview screen and perform manipulations such as zoom, PAN and rotate for easier image diagnosis. Switch to the viewer screen to review all patient studies and access extensive toolbar.



PLAN your studies



POSITION your patient



PREVIEW your image



A COMPLETE VETERINARY SPECIFIC DIGITAL IMAGING SYSTEM

SPECIFICATIONS

DelWorks DR Detector

- TOSHIBA FDX4343R, Cesium Iodide (CsI) Scintillator with Amorphous Silicon (a-Si) Photodiode
- 17in x 17in detector size (512mm x 495mm x 43mm)
- 16.9in x 17.3in active area
- 143µm pixel size
- 9 megapixel (3,008 x 3,072 pixels)
- 3.5 lp/mm
- 6 second cycle time (shot to shot)
- MTF (2.0 lp/mm, 70 kVp, 1x1) 30% typ.
- DQE (0), Quantum-Limited >70%
- . Less than one (1) second to image preview
- Operating temperature +50° F to +95° F
- 14 bit A/D conversion
- 16 bit digital output (1000 Base-T)
- External X-ray synchronization control using included external synchronization device

Technologist Workstation Hardware

- Dell Optiplex 7010, small form factor computer with Intel Core i5, 3.0GHz, 4GB RAM, 500GB hard drive, weight 15lb (6.80kg)
- 19in touch screen, flat-panel display, weight 18.74lb (8.5kg). Display rotates 90° landscape/vertical. Active screen diagonal 19in (483mm):
- width 14.80in (376mm); height 11.85in (301mm) Pixel pitch 0.294mm, 16.8 million colors
- Keyboard and mouse

Processing Software

- Extensive measuring functions
- CD/DVD burning for patient images
- DICOM Store, Print, and MWL
- Generator integration (on select models)
- Advanced Image Processing features
- Unique Animal Positioning Guide
- Grid suppression software
- Extensive tool library with measuring,
- layout adjustments, export options and more

16.46in (418mm

(70mm)

Patient CD creation

3 654



3 5in

(9cm)



DEL MEDICAL

Toll Free: (800) 800-6006 (847) 288-7011 Fax: Email: info@delmedical.com www.delmedical.com 241 Covington Drive, Bloomingdale, IL 60108

Universal

Universal branded products are manufactured and distributed by UMG/Del Medical. Del Medical is an ISO 13485 & ISO 9001 Certified facility. Specifications are subject to change without notice.

Toll Free: (800) 261-9808

(8.9cm) (8.9cm)

Fax: (914) 835-6111

Email: info@umgxray.com

3 5in

www.umgxray.com

28 Calvert Street, Harrison, NY 10528

DVChoice Integrated Table/Tube Stand

- · Welded table construction with 400lb load capacity and urine trap
- Integrated tube stand with electric locks travels full length of table
- · Angulating tube arm, angulation dial and operator handle
- Electric table locks on 2-way float or 4-way float table top
- Foot activated exposure switch
- 17in x 17in grid cabinet with tube stand interlock for FDX4343R detector
- Grid: 10:1 ratio, 178 line, 34-44in focal range
 - Certified manual collimator with built-in spare bulb and retractable tape measure
- X-ray tube: 200,000 heat unit, 0.6-1.5mm focal spots with 90° horn angle
- · One pair 15ft high voltage cables, terminated both ends

DVChoice System Options

· Available with stationary, 2-way float or 4-way float table top

Generator

 CMDR 500mA, 40kW, high frequency, digital ready, three-phase with anatomical programming

Generator Options

- Pedestal stand
- 32kW generator options are also available

Side View **DVChoice all models**



3.5in

(8.9cm)

68.5in

(174cm)

1 17 (79cm)



G

. 1ω

#9000-DVChoice

3 5in