

### Technology

Through our commitment to quality and innovation, DEL MEDICAL continues to extend limits and maintain a leading role in the development of Healthcare products that better the needs and wants of global communities. It is our goal to exceed expectations by collaborating with our customers and partners to focus on advanced technologies and inspire enthusiasm with high-value-added products.

### Experience

When you make the choice to partner with DEL MEDICAL, you join one of the most experienced medical imaging teams in the industry. For over 85 years we have been dedicated to providing customers with the highest quality products and services. We strive to improve work flow and create further efficiencies while enhancing your bottom line.

### Support

DEL MEDICAL equipment products are serviced through local, trained professionals around the globe. In addition, we invest a lot of technical and human power to offer expert support assistance from a dedicated team. Located throughout North America the DEL MEDICAL Customer Assistance team is among the most knowledgeable in the industry and are ready to provide fast, courteous troubleshooting assistance.



### www.delmedical.com



**New York** 28 Calvert Street









#### 241 Covington Drive Harrison, NY 10528 Bloomingdale, IL 60108 (Tel) 800.800.6006 (Tel) 800.261.9808 (Fax) 914.835.6111 (Fax) 847.288.7011

Universal FINTS DIGITAL RADIOGRAPHIC FLOOR MOUNTED SYSTEM





## Affordable quality and reliability

The FMTS Digital Radiographic System is powered by our DelWorks T-Series DR, combining the in-room flexibility and portability of a wired 14 in x 17 in (35 cm x 43 cm) Cesium-lodide digital imaging detector with the powerful diagnostic usability of the radiographic components.

Add an additional detector such as a 17 in x 17 in (43cm x 43cm) fixed-mounted detector into the Wall Stand and further increase the efficiency and productivity of your new Del Medical Digital Radiographic Imaging Suite.



## FMTS Tube Stand Precise Accuracy

- Exceptional floor mounted movement in low-ceiling environments
- Functionally designed tube head hand grip for quick, proper positioning
- All locks release button for smooth, multi-directional movement
- Column rotation for cross table lateral exams

## EV650 Table Versatile Comfort

- 650 lb. (294 kg) patient weight capacity
- Heavy-duty electric table locks with foot treadle release
- Lowers to 22.25 in (56.5 cm) for comfortable patient transfer
- Various table length and width options
- Stain resistant fiber resin top

## VS200 Wall Stand Upright Convenience

- Effortless, vertical movement for fast patient positioning
- Patient-side and optional overhead handgrips for patient stability and repeatable positioning
- Receptor lowers to the floor offering a full range of upright weight bearing procedures

# Full generator control and Cesium detector options

The DelWorks T-Series DR System provides a user friendly single-point of operation for your system compatible generator, along with integration to your patient information system for a seamless workflow experience.

### **T14C Wired Detector**

The T14C wired, 14 in x 17 in (35cm x 43cm) detector is ultra-light and fits perfectly inside the VS200 wall stand and EV650 table bucky without any modifications. The in-room portability of the T14C allows the versatility needed to optimize workflow without the need for batteries.

## **T17C Fixed Detector**

The T17C fixed, 17 in x 17 in (43 cm x 43 cm), large format flat-panel detector is designed to be mounted inside the EV650 table and/or VS200 wall stand. It delivers high resolution, low noise and superior image quality along with the benefit of lower radiation dose to the patient.





